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CLSI. M100-S26 performance standards for antimicrobial susceptibility testing; Twenty-sixth informational supplement; 2016.

CLSI. M45-S3 Methods of antimicrobial dilution and disk susceptibility testing of infrequently isolated or fastidious bacteria; Third edition; 2016.

Humphries RM, Linscott AJ. Laboratory diagnosis of bacterial gastroenteritis. Clin Microbiol Rev. 2015 Jan;28(1):3-31. (formerly Cumitech)

Acheson DWK. Differential diagnosis of microbial foodborne disease. www.uptodate.com, update 08-05-2017.

Barshak MB. Antimicrobial approach to intra-abdominal infections in adults. www.uptodate.com, update 14-03-2016.

Benagli C, Demarta A, et al. A rapid MALDI-TOF MS identification database at genospecies level for clinical and environmental *Aeromonas* strains. PLoS One. 2012;7(10):e48441.

Chen PL, Tsai PJ, et al. *Aeromonas* stool isolates from individuals with or without diarrhea in southern Taiwan: Predominance of *Aeromonas veronii*. J Microbiol Immunol Infect. 2015 Dec;48(6):618-24.

Darby A, Lertpiriyapong K, et al. Cytotoxic and pathogenic properties of *Klebsiella oxytoca* isolated from laboratory animals. PLoS One. 2014 Jul 24;9(7):e100542.

Gelfand MS. Clinical manifestations and diagnosis of *Listeria monocytogenes* infection. www.uptodate.com, update 16-06-2016.

Harris JB, LaRocque RC, et al. Cholera. Lancet. 2012 Jun 30;379(9835):2466-76.

Hirai Y, Asahata-Tago S, et al. *Edwardsiella tarda* bacteremia. A rare but fatal water- and foodborne infection: Review of the literature and clinical cases from a single centre. Can J Infect Dis Med Microbiol. 2015 Nov-Dec;26(6):313-8.

https://nrchm.wiv-isp.be/nl/ref_centra_lab/shiga_toxine_verotoxine/default.aspx

<https://www.cdc.gov/narms/pdf/2014-Annual-Report-narms-508c.pdf>

<http://www.eucast.org/>

Janda JM, Abbott SL. Infections associated with the genus *Edwardsiella*: the role of *Edwardsiella tarda* in human disease. *Clin Infect Dis*. 1993 Oct;17(4):742-8.

Janda JM, Abbott SL, et al. *Plesiomonas shigelloides* Revisited. *Clin Microbiol Rev*. 2016 Apr;29(2):349-74.

Kaiser L, Surawicz CM. Infectious causes of chronic diarrhoea. *Best Pract Res Clin Gastroenterol*. 2012 Oct;26(5):563-71.

Lamy B, Kodjo A, et al. Identification of *Aeromonas* isolates by matrix-assisted laser desorption ionization time-of-flight mass spectrometry. *Diagn Microbiol Infect Dis*. 2011 Sep;71(1):1-5.

Lamy B, Kodjo A, et al. Prospective nationwide study of *Aeromonas* infections in France. *J Clin Microbiol*. 2009 Apr;47(4):1234-7.

Lamy B, Laurent F, et al. Accuracy of 6 commercial systems for identifying clinical *Aeromonas* isolates. *Diagn Microbiol Infect Dis*. 2010 May;67(1):9-14.

Lobatón T, Hoffman I, et al. *Aeromonas* species: an opportunistic enteropathogen in patients with inflammatory bowel diseases? A single center cohort study. *Inflamm Bowel Dis*. 2015 Jan;21(1):71-8.

Man SM. The clinical importance of emerging *Campylobacter* species. *Nat Rev Gastroenterol Hepatol*. 2011 Oct 25;8(12):669-85.

Manual of Clinical Microbiology, 11th Ed. 2015, Versalovic J (Eds.), ASM Press Washington DC.

Noonan L, Freeman J. *Bacillus cereus* and other non-anthraxis *Bacillus* species. www.uptodate.com, update 09-02-2017.

Ooi ST1, Lorber B. Gastroenteritis due to *Listeria monocytogenes*. *Clin Infect Dis*. 2005 May 1;40(9):1327-32.

Ooka T, Seto K, et al. Clinical significance of *Escherichia albertii*. *Emerg Infect Dis*. 2012 Mar;18(3):488-92.

Pegram PS and Stone SM. Botulism. www.uptodate.com, update 07-05-2016.

Perry JD. A Decade of Development of Chromogenic Culture Media for Clinical Microbiology in an Era of Molecular Diagnostics. *Clin Microbiol Rev*. 2017 Apr;30(2):449-479.

Pfeiffer ML, DuPont HL, et al. The patient presenting with acute dysentery--a systematic review. *J Infect*. 2012 Apr;64(4):374-86.

Raszl SM, Froelich BA, et al. *Vibrio parahaemolyticus* and *Vibrio vulnificus* in South America: water, seafood and human infections. *J Appl Microbiol*. 2016 Nov;121(5):1201-1222.

Schlech WF 3rd, Schlech WF 4th, et al. Does sporadic *Listeria* gastroenteritis exist? A 2-year population-based survey in Nova Scotia, Canada. *Clin Infect Dis*. 2005 Sep 15;41(6):778-84.

Shah MM, Odoyo E, et al. First Report of a Foodborne *Providencia alcalifaciens* Outbreak in Kenya. *Am J Trop Med Hyg*. 2015 Sep;93(3):497-500.

Shima A, Hinenoya A, et al. Prevalence of *Providencia* Strains among Patients with Diarrhea and in Retail Meats in Thailand. *Jpn J Infect Dis*. 2016 Jul 22;69(4):323-5.

Shin HB, Yoon J, et al. Comparison of MALDI-TOF MS, housekeeping gene sequencing, and 16S rRNA gene sequencing for identification of *Aeromonas* clinical isolates. *Yonsei Med J*. 2015 Mar;56(2):550-5.

Zollner-Schwetz I, Högenauer C, et al. Role of *Klebsiella oxytoca* in antibiotic-associated diarrhea. *Clin Infect Dis*. 2008 Nov 1;47(9):e74-8.